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This warranty certificate entitles either "free-of-charge" replacement or repair of the components recognized faulty by the manufacturer, provided the instrument is returned "freight-pre-paid" to the distributor.

The following are excluded:

- **ALL GLASS-WARE**
- Damages caused during the shipment
- Faulty power connections
- Accidental unforeseen faults
- Damages caused by improper use and/or not in conformity with the enclosed instructions
- Product liability is revoked if anyone tampers with the instrument or its additional parts anyway

FOR RESELLER
ONLY

Purchaser's name
Purchaser's address
Zip City Country.....
Date of purchase..... Seller's name.....

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NOTE

**Read the instructions carefully before operating the Balance.
This is an accuracy instrument and must be handled with care.**

INDEX

WARRANTY

WARNING

DECLARATION OF CONFORMITY

1	WARNING	Pag. 23
2	TECHNICAL CHARACTERISTICS.....	Pag. 24
3	LOCATION	Pag. 24
4	INSTALLATION	Pag. 25
5	OPERATING CONTROLS OF THE BALANCE.....	Pag. 26
6	OPERATING INSTRUCTIONS.....	Pag. 27
7	KEYBOARD FUNCTION.....	Pag. 28
8	TARE	Pag. 29
9	MAIN MENU AND OPERATING FUNCTIONS.....	Pag. 29
9.1	UNIT2	Pag. 30
9.1.1	STANDARD BALANCE FUNCTION	Pag. 30
9.1.2	PIECE COUNTING FUNCTION.....	Pag. 30
9.1.3	PERCENTAGE FUNCTION.....	Pag. 31
9.2	SET UP	Pag. 32
9.2.1	SERIAL	Pag. 33
9.2.1.1	BAUD	Pag. 33
9.2.1.2	BYTE FOR	Pag. 33
9.2.1.3	PROTOC	Pag. 33
9.2.1.4	TRANS OPT	Pag. 33
9.2.2	CALIBR.....	Pag. 34
9.2.2.1	CAL OFF.....	Pag. 34
9.2.2.2	CAL ON	Pag. 34
9.2.2.3	CAL AUTO	Pag. 34
9.2.3	A-ZERO	Pag. 34
9.2.3.1	ZERO OFF	Pag. 34
9.2.3.2	ZERO ON	Pag. 34
9.2.4	FILTER	Pag. 34
9.2.4.1	SLO	Pag. 34
9.2.4.2	AVG	Pag. 34
9.2.4.3	FAST	Pag. 34
9.2.5	SET-DEF	Pag. 34
9.2.5.1	END	Pag. 34
9.2.5.2	STORE	Pag. 34
9.3	CALIBRATION	Pag. 35
9.3.1	EXTERNAL CALIBRATION	Pag. 35
9.3.2	INTERNAL CALIBRATION	Pag. 35
9.3.3	TIMED AUTOMATIC CALIBRATION	Pag. 35
10	QUICK RESTORE OF DEFAULT PARAMETERS.....	Pag. 36
11	SERIAL OUTPUT.....	Pag. 37
11.1	SERIAL OUTPUT CONNECTIONS	Pag. 37
11.2	DATA TRASMISSION FORMATS.....	Pag. 37
11.2.1	CRYSTAL protocol	Pag. 37
11.2.2	EURO protocol	Pag. 38
11.2.3	PRINTER protocol	Pag. 38
11.3	SERIAL COMMANDS.....	Pag. 39

1 WARNING

The power supply pack of the CENT balance complies with Electrical Protection Class II.
Before connecting the power, check that the operating voltage, indicated on the Power supply pack, matches the line voltage (100 / 240 Vac).

CENT BALANCE MAY ONLY BE USED WITH THE POWER SUPPLY PACK.

THE BALANCE MUST BE OPERATED ONLY IN DRY, CLEAN AND IF POSSIBLE VIBRATION-FREE ENVIRONMENT.

THE BALANCE MUST NOT BE OPERATED IN AREAS WHERE THERE ARE EXPLOSION RISKS.

NOTE:

In case of improper setting of the balance parameters during setup, the balance may power-up in an incorrect configuration.

The user can restore the Factory configuration to avoid the malfunction as follows:

Disconnect the balance supply.

Keep the On/OFF key pressed

Reconnect the supply

Release the On/OFF key when the display shows 'Set Factory'

Now the Factory default settings are restored.

2 TECHNICAL CHARACTERISTICS

	MODELS					
	CENT 2000	CENT 2000 CAL	CENT 4000	CENT 4000 CAL	CENT 6000	CENT 6000 CAL
Capacity and readability	2200 g / 0,01 g		4200 g / 0,01 g		6300 g / 0,1 g	
Linearity	± 0,02 g				± 0,2 g	
External calibration weight	500 g or multiplies		1000 g or multiplies			
Average stabilization time: fast / slow (seconds)	2 / 4					
Operating temperature	18 - 28 °C (10 - 30 °C at reduced accuracy)					
Power Voltage	100 / 240 Vac -15/+10% (range from 85 to 264 V) 50/60Hz					
Power consumption	200 mA					
Data Output	EIA standard RS232					
Balance dimensions	216 x 380 x 125 (b x p x h)					
Pan type	Ring			Ring		
Pan dimension (mm)	160 x 160			160 x 160		
Balance weight	<u>5,3 Kg without internal calibration mechanism</u> <u>6,1 Kg with internal calibration mechanism</u>					
Internal calibration	NO	YES	NO	YES	NO	YES

3 LOCATION

Choose an appropriate location taking into account these criteria:

- ◆ Vibrations and air drafts in the ambient must be avoided.
- ◆ The supporting stand must be rigid.
- ◆ The main line must be stable. It is not advisable to use the same electric line of the fluorescent lamps or other electric apparatus that can be switched on and off randomly.
- ◆ The balance must be kept far from magnetic fields.
- ◆ Environmental conditions with wide temperature changes must be avoided. Ambient temperature above 28 °C or below 18 °C may affect the operation of the balance and its accuracy.


In very dry environments electrostatic charges may build up, especially if glass or porcelain objects are used

4 INSTALLATION


Instructions to remove the locking device of the balance

Unpack the balance and all its accessories.

- ◆ Unscrew the knob on the locking bracket as indicated by the arrow (anti clockwise).
Retain the locking bracket in case the balance has to be returned to the supplier. This bracket can be refastened to the balance by reolacing the knob.
- ◆ Screw the pan with the central screw
- ◆ Level the balance using the spirit level on the rear and rotating the levelling feet under the front part of the basement (Fig 2).
- ◆ Insert the power supply cable located on the rear panel (Fig 2).
- ◆ Insert the power supply in the main line plug checking the supply voltage printed on the power supply case (100 / 240 Vac) - (Fig 2).

Rear panel details - FIG.2

5 OPERATING CONTROLS OF THE BALANCE

FIG. 3

Weighing pan
Display
Print Key
Mode Key
Tare Key
ON-OFF Key

6 OPERATING INSTRUCTIONS

Wait for about 15 minutes after power up before performing any weighing operation.
It's advisable to leave the balance permanently connected to the Ac line.

Press the ON/OFF key more than one second. The LCD backlight will be turned off and on the display will appear the message "OFF".

7 KEYBOARD FUNCTION

The Cent's keyboard has four key labeled: PRINT, MODE, TARE and ON-OFF.

Every key has at least one function (called short function) while key MODE and ON/OFF also activate a second function (called long function)

To activate short function presses the key and release it before one second.

To get long function press the key and keep pressed for more than one-second .

In the following table are shown the all keyboard functions.

Normal Use

KEY NAME	KEY ACTIVATED FOR LESS THAN ONE SECOND (short function)	KEY ACTIVATED FOR MORE THAN ONE SECOND (long function)
On/OFF	Display illumination ON – OFF	Balance ON - OFF
Mode	See additional details	See additional details
Tare	Tare function	Tare function
Print	Print function	No function

Main menu and configuration details

KEY NAME	KEY ACTIVATED FOR LESS THAN ONE SECOND (short function)	KEY ACTIVATED FOR MORE THAN ONE SECOND (long function)
On/OFF	Confirm selection	No function
Mode	Scroll available items	No function
Tara	Exit current selection	No function
Print	No function	No function

Keyboard details - FIG. 4

8 TARE

After a container is placed on the pan, press the TARE key.

The LCD display will show zero value.

The balance has a tare capacity up to its measuring range.

Place the sample to be weighed in the container.

The display shows the sample weight.

It's possible to use the tare function after adding several samples up to the maximum balance capacity indicated by a bargraph on the left side of the LCD display (FIG.4).

9 MAIN MENU AND OPERATING FUNCTIONS

The balance can be operated in the following configurations:

Balance with gram, carat, ounce, pound measuring unit

Piece counting balance

Percentage computing balance

The selection is performed as follows:


Hold MODE key pressed until the display shows "UNIT 2", release the key. Pressing key MODE (short time) it is possible to select the functions "UNIT 2", "SETUP", "CAL E", wrapping around at the end of the list.

The ON/OFF key activates the selected submenu, while the TARE key exits the main menu.

When the submenu is entered, the keys operate as above: pressing ON/OFF the selection is confirmed, while pressing TARE the selection is left unchanged.

NOTE: If no key is activated for 20 seconds, the menu level is incremented up to the end, without altering any parameter and the balance resumes its normal operation.

Main Menu flow

UNIT 2

Submenu defining the additional operating functions of the balance.

SETUP

Function used to configure the balance


CAL E

Activation of the external calibration (manual calibration with user supplied weight)

CAL Int

activation of the internal calibration (automatic calibration with built-in weight)

9.1 UNIT 2

With this function the balance can be programmed to operate with different measuring units. CENT can weight also in ounce (OZ), pounds (LB) and carats (CT). The balance can also be programmed for operation as:

- ◆ STANDARD BALANCE
- ◆ PIECE-COUNTING BALANCE
- ◆ PERCENTAGE BALANCE

Select the desired unit or operating mode, using the MODE key.
Confirm the unit/function selected with the ON-OFF key.
Now the balance is operating with the new unit or function.

9.1.1 STANDARD BALANCE FUNCTION

When not operating as piece counting or percentage balances the keys act in normal mode. The key MODE (activated for less than 1 second) switches between gram and the second unit. If a measuring unit different from gram (main operating unit) is selected, the number of decimal places is changed and the resolution of the balance (total number of divisions) is reduced, so that the last digit always exhibits a 1-count increment. When a second unit is selected the balance automatically stores the information so that, even if the power is removed, the selection is retained and, when the balance is switched on, the selected unit will be active.

9.1.2 PIECE COUNTING FUNCTION

The function is selected as above with the menu item PCS. Press the key mode more than one second. Display will show 5 PCS blinking: 5 PCS is the first sample number for acquisition. If another sample is desired, then press the key mode repeatedly. The display will show 10 PCS, 20 PCS, 50 PCS and finally 100 PCS. Then, the PCS menu rolls over and 5 PCS message is shown again.

When the reference sample is on the pan and the stability signal light is on, push the ON-OFF key to confirm sample acquisition. Automatically the new reference value is acquired and balance starts to display the piece count.

Note:

If the weight of the sample is too small, acquisition is interrupted and the message "PCS ERROR" appears until the TARE key is activated.

9.1.3 PERCENTAGE FUNCTION

The function is selected as above with the menu item %.

A new reference may be acquired pressing the MODE key for more than 1 second.

The display will show – 100 – blinking.

Now the reference quantity may be weighed, pressing the ON-OFF key, when the stable symbol lights on the display, the new reference value is acquired and the balance starts to display the weight as a percentage of the reference acquired.

Notes:

If the weight of the sample is too small, acquisition is inhibited and the message “PERC ERROR” appears until the TARE key is activated.

The balance can measure percentages up to 500 % of the reference if the full scale is not exceeded.

Example for percentage calculation:

If a reference weight of 180 g is acquired the display shows 100.00 %. Now, if a weight of 35 g is placed on the pan the display shows 19,44 %. If the weight is increased up to 95 g, the display will show 52,77 %.


Now, if the weight is increased more than 900 g the display will show the message OVER PERC.

This message means that whole percentage function capacity (500 %) has been exceeded.

Pressing the key MODE, you can switch from main unit to percentage unit to know the weight value loaded on the pan.

To go out from OVER PERC condition, remove the weight from the pan.

9.2 SET-UP

9.2.1 SERIAL

This menu controls the serial output control parameters.

9.2.1.1 BAUD

It defines the transmission baud rate of the serial output from 1200 to 9600 baud.

9.2.1.2 BYTE FOR

It defines the format of the transmitted byte and may be set as:

DATA FORMATS	DATA BIT	PARITY CHECK TYPE	BIT STOP
7 E 1	7	EVEN	1
7 O 1	7	ODD	1
7 N 2	7	NO PARITY	2
7 E 2	7	EVEN	2
7 O 2	7	ODD	2
8 N 1	8	NO PARITY	1
8 E 1	8	EVEN	1
8 O 1	8	ODD	1
8 N 2	8	NO PARITY	2

9.2.1.3 PROTOC

It is the protocol of the serial data output and may be:

EUROPE
CRYSTAL
PRINTER

For further details refer to DATA TRANSMISSION section

9.2.1.4 TRANS OPT

9.2.1.4.1 REQUEST

Data is transmitted on-demand data transmission

9.2.1.4.2 CONTINUE

Data is transmitted every balance conversion.

It may be used to drive peripherals than can't issue a print command to the balance.

9.2.1.4.3 DATA STABLE

Data is transmitted if and only if the weight on the pan is stable.

These functions can be stored by STORE function (see 9.2.5)

9.2.2 CALIBR

The CENT balance may be calibrated with an external mass.

9.2.2.1 CAL OFF → Both calibration are disabled

9.2.2.2 CAL ON → Both calibration are enabled

9.2.2.3 CAL AUTO → With this selection the calibration is enabled and a timer activates the internal calibration every 3 hour. The automatic calibration is delayed if, during last 10 minutes, the balance has been operated or serial comand has been issued.

9.2.3 A-ZERO

If the Auto-Zero function is active and when no weights are on the pan, the automatic zero tracking correction is active: the balance remains at a net value equal to zero. The function is useful to cancel small drift of the zero reading.

9.2.3.1 ZERO OFF

The zero tracking function is inactive.

9.2.3.2 ZERO ON

The function is active and small changes in zero reading are corrected.

9.2.4 FILTER

The setting of the filter menu controls the dynamic response of the balance. Selections are.

9.2.4.1 SLO slow response applicable to disturbed environments (high accuracy weighing)

9.2.4.2 AVG average value for both high accuracy and

9.2.4.3 FAST fast response applicable to quickly changing weights (dosing)

9.2.5 SET DEF

All settings, changed by the setup menu are valid as soon as they are changed.

However if the balance is powered off they are lost and when powered up the balance will use different parameters.

To maintain a custom configuration that will be active at power-on the function SET DEF is implemented.

9.2.5.1 STORE → current settings are stored as user default

9.2.5.2 END → user default settings are left unchanged

9.3 CALIBRATION

9.3.1 EXTERNAL CALIBRATION

Once the external calibration function is started the display shows "LOAD" and the reference weight acquisition is active: a reference weight of E2 class is required for calibration, with value function of the range, the value of the weight is automatically detected and checked by the balance.

When the value of the reference weight is acquired, the display shows the "UNLOAD" message: at this time the weight must be removed from the pan and the balance enters the zero acquisition phase.

At the end of the cycle the message "CALIBRATED" appears if the calibration is valid and the new calibration factor is stored, if anything fails, the message "CAL ERROR" appears and the calibration factor of the balance remains unchanged.

Notes:

see section 2 for the value of the calibration mass

Note: The calibration process not starts if a weight of more than twenty grams are on the pan at calibration start time.

9.3.2 INTERNAL CALIBRATION

If the balance is equipped with the internal calibration mass the message "CAL INT" appears as a choice of the main menu and the internal automatic calibration is available.

Once this function is activated, the automatic calibration process begins.

The display shows "CAL AUTO" during the calibration process while the functions of reference and zero weight acquisition are performed automatically by the balance.

During the calibration the process is monitored for proper execution.

At the end of the cycle the message "CALIBRATED" or "CAL ERROR" appear as a result of the calibration process as above.

9.3.3 TIMED AUTOMATIC CALIBRATION

If the balance is equipped with the internal calibration mass, the timed calibration periodically starts an internal automatic calibration.

In order to operate, 3 hours must be elapsed from the last calibration, the balance must be inactive for at least 10 minute before the calibration time, no serial command was send to the balance.

The display shows "CAL AUTO" during the calibration process while the functions of reference and zero weight are performed automatically by the balance.

At the end of the cycle the message " CALIBRATED " or "CAL ERROR" appear as a result of the calibration process as above.

In this way the operator is able to determine the status of the last timed calibration.

To restore the normal operation the key TARE must be activated.

10 QUICK RESTORE OF DEFAULT PARAMETERS

If configuration parameters are modified, so that a wrong setup is stored, a quick restore function of the factory parameters is available.

Remove the power supply plug from the rear of the balance

Press and keep pressed the ON/OFF key

Reinsert the power supply plug

Wait until the display shows "CENT"

Release the ON/OFF key

The display will show "SET FACTORY"

Factory default parameters are restored

Note:

When the balance configuration has been modified and its operation is correct, the use of the SET DEF command of the setup menu will permanently store the new configuration parameters in memory so that they are not lost on power down.

11 SERIAL OUTPUT

CENT serial output is compatible with EIA RS232 standard.

Every answer from balance is sent every 100 ms.

11.1 SERIAL OUTPUT CONNECTION

BALANCE	DB9 CONNECTOR (9 pole)	DB25 CONNECTOR (25 pole)
RS 232 Out	Pin 2	Pin 3
RS 232 In	Pin 3	Pin 2
GROUND	Pin 5	Pin 7

11.2 DATA TRANSMISSION FORMAT

The format of the data transmitted by the balance is different as a function of the selected protocol by the menu option 9.2.1.3.

EURO and CRYSTAL protocols are suitable for computer interfacing; PRINTER protocol may be used to drive external printers.

11.2.1 CRYSTAL protocol

The answer to any type of command is formatted with a data field, a unit field and a status field as follows:

a measure field
a unit field
a status field

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	B	F1	F2	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19

The measure field (1..10) is the reading on the display.
It is right justified with decimal point and measure sign.

Unit field shows the current unit weight

Status field show the current balance status

B is the blank character.

Field N°18 and 19 are:

CR carriage returns ASCII code.

LF line feed ASCII code.

Status information of the F1 character

D	Transmitted data is valid
O	Balance in over range
U	Balance in under range
T	Tare
C	Calibration process in progress
Z	Zero acquisition
I	Init test

Status information of the F2 character

S	Measure is stable
I	Measure is not stable
E	Error
A	Automatic calibration (when F1 = C)
L	Load weight(manual calibration) (when F1 = C)
U	UnLoad weight (manual calibration) (when F1 = C)
D	Calibration done (when F1 = C)
B	Calibration busy(already started)
O	Calibration OFF
P	Percentage

11.2.2 EURO protocol

The answer to any type of command is formatted as follow:

D0	D1	D2	D3	D4	D5	D6	D7	D8	CR	LF
1	2	3	4	5	6	7	8	9	10	11

The measure field (1..9) is the display reading.
It is right justified with decimal point and the measure sign.

Field N°10 and 11 are:

CR carriage returns ASCII code.
LF line feed ASCII code.

11.2.3 PRINTER protocol

Is activated by the "PRINT" pushbutton. The output string uses 10 characters for the measure including sign and decimal point, the value is right justified and filled with spaces from the left. The format of the output is:

D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	B	U	U	U	CR	LF
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16

The measure field (1..10) is the display reading.
It is right justified and has the decimal point and the measure sign.

Unit field shows the current unit weight

B is the blank character.

Field N°15 and 16 are:

CR carriage returns ASCII code.

LF line feed ASCII code.

Note

Data will be printed if and only if the weight on pan is stable.

11.3 SERIAL COMMANDS

The balance responds to the following commands only if a CRYSTAL or EURO protocol is selected:

SERIAL REMOTE CONTROLS TABLE

COMMAND	FUNCTION
B	Send data on display
C	Start calibration
I	Start continuous transmission of data on display
F	Stop continuous transmission of data on display
T	Tare command

Note 1:

Every remote control command must end with <CR>.

Line-feed characters are ignored.